

Rider 50a Endoscope Washer-disinfector

2MED RIDER50A



Product Characteristics

The most speedy and wash disinfection Touch-screen so that the program parameters can be displayed and controlled Auto door with foot switch Automatic visual glass door adopted for observing the operation inside the tank Totally enclosed to avoid sanitizer odor leakage Automatic adjuvant injection: only need to link the adjuvant container with the tank through a connector Work with various kinds of sanitizer: glutaraldehyde, peracetic acid, phthalic dicarboxaldehyde, chlorine dioxide and acidic electrolyzed oxidizing water.

Includes leaking test, blocking test, wash, enzyme wash, rinse, disinfection, final rinse, alcohol drying and air drying Data management and printing for disinfection information Overall dimensions: 74x65x94 Weight 100Kg Water Consumption/cycle: 25L Dynamic drainage

Acail Angola Polo Industrial de Viana EN Luanda - Catete Km 23 ANGOLA

922 327 152 acailangola@gmail.com www.acailangola.com